# lumie<sup>®</sup> zip<sup>™</sup>

**OPERATING INSTRUCTIONS  
NOTICE D'UTILISATION  
BEDIENUNGSANLEITUNG  
BRUKSANVISNING  
GEBRUIKSAANWIJZING**

English

## OPERATING INSTRUCTIONS

### Zip

---

- 1 Precautions
- 2 Safety
- 3 User guide
- 4 When should I use it?
- 5 What time of day?
- 6 How long?

English

## INTRODUCTION

Using Zip will provide you with the bright light your body needs during the dark winter months. It can help to boost your energy levels, put you in a better mood and make you feel more awake. Light therapy is safe and is proven to help alleviate SAD (Seasonal Affective Disorder) and the winter blues\*.

There is growing evidence that our circadian rhythms are particularly sensitive to light from the short wavelength part of the spectrum. Traditional, large lightboxes use fluorescent tubes and offer high levels of brightness (usually measured in lux) across all wavelengths of the light spectrum. LEDs lend themselves to smaller, more portable light therapy designs with a corresponding reduction in brightness. To help make up for the lack of lux, the LEDs in your Zip have been specially selected to emit lots more light at the specific short wavelengths that help to regulate melatonin, making it a handy, as well as effective, light therapy option.

Please read these instructions carefully to get the most out of Zip.

\*Am J Psychiatry. 2005 Apr; 162(4):656-62

English

## 1 Precautions

Research shows that light therapy is very safe. There are no problems with ultraviolet (UV) light as Zip emits only a tiny amount, far less than normal daylight. However, please check with your doctor first if any of the following apply:

- You have a major eye problem
- You have been advised to avoid bright light
- You are taking medicine where the instructions advise you to avoid bright light
- You have, or have had, deep depression.\*

\*This is especially important if you are also using medication, since Zip is designed to help improve your mood

If you have any problems whilst using Zip, or if you've been using the light for a week but haven't noticed any benefit, please call the Lumie Careline on 01954 780500 or speak to your doctor.

**IMPORTANT:** the guidelines in this leaflet are not medical directions for treating any condition. If your doctor has suggested using light therapy, follow their instructions and discuss any effects with them.


## 2 Safety

- Keep the unit away from heat, water and damp. For indoor use only.
- Children using the light, or playing near it, should be supervised.
- Use only the adapter supplied.
- Do not attempt to service any part of the unit or replace the battery.
- Do not dispose of the battery in a fire or incinerator.

English

### 3 User guide

Zip can be powered from the mains or the internal rechargeable battery. The battery takes around two hours to fully charge and then provides around two hours use. To use the mains supply or to charge the battery, plug in the adapter and switch on at the mains.

- Press the side button to turn on the battery charge indicator: this flashes when charging.
- Press again to turn on the light and activate the timer.
- Press the button as necessary to set the timer: display shows 10-minute segments up to 60 minutes.
- Timer segments disappear one by one and at the end of your treatment the light will automatically fade down and turn off.
- To turn off the light at any time, press and hold the side button. On release, the light will fade down and turn off.
- The battery charge indicator will remain on for around 15 minutes.

## English

You shouldn't stare directly at the light but it must reach your eyes to have any effect. Position Zip at arm's length (about 50cm) and use the stand to angle the light towards your face. You can use Zip whilst you're on the computer, when you're watching TV, at mealtimes, etc. Dark glasses, tinted lenses or falling asleep will all limit the effects of light therapy.

You should notice positive benefits after using Zip for 3 or 4 days in a row.

## 4 When should I use it?

It's best to tackle SAD or the winter blues as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable, start using your Zip. If you learn to use the light early on then you should never dread the winter months again.

Most people can skip a day here and there, especially when the weather is bright. As spring comes around you can probably start to cut back your light therapy; your body will soon tell you if it's too soon and if you notice your symptoms returning simply use your Zip again for a few days.

English

## 5 What time of day?

It depends on your sleep cycle but once you've found a pattern that's right for you, try to use it around the same time every day.

### **I find it difficult to get going in the morning**

Switch on Zip as soon as you can. If you don't have time in the morning, then you should aim for a top-up or full session in the afternoon. You could take it to work if you're at a desk most of the day.

### **I feel sleepy or fall asleep too early**

Use Zip in the late afternoon/early evening to perk you up. It's best to avoid bright light too close to bedtime though - within about three hours - as you may find it difficult to get to sleep.

## 6 How long?

Everybody's different but, if you use Zip at about 50cm, as suggested, around 30 minutes will probably be enough. If you're more comfortable with it further away, then less light will reach your eyes and you'll need to use it for longer. If you find your eyes ache slightly, try moving it further away or switch off for a while.

English

## Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you are returning the unit for service, it must be in its complete original packaging including internal pieces. Make sure that the mains plug is prevented from moving in transit and damaging the unit. If the unit is received damaged or without original packaging, then we may have to charge to put it right. Please call us for advice before sending the unit back.

## Technical specification

Subject to change without notice.

Lumie Zip

Input 100-240Vac 0.4A

Output 12Vdc 1.00A max

Operating temperature 0°C-35°C

Class IIa Medical Device

English

## Contact us

Lumie is a trademark of Outside In, Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits. If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on 01954 780500 or email [info@lumie.com](mailto:info@lumie.com).

Lumie

3 The Links, Trafalgar Way

Bar Hill

Cambridge

CB23 8UD

UK

Lumie is a trademark of:

Outside In (Cambridge) Limited incorporated in England and Wales.

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© Lumie 2008
Outside In (Cambridge) Limited met rechtspersoonlijkheid in Engeland en Wales.
Outside In (Cambridge) Ltd