lumie®

halo™

 ENG | FRA | ESP

Safety

WARNING - To reduce the risk of fire, electric shock or injury:

Read instructions for use. 

If the unit is damaged, contact Lumie before use.

Keep away from water and damp. 

For indoor use only. 

The unit gets warm when the light is on, so make sure that it is not covered.

Use only the mains power adaptor supplied.

Keep the mains power adaptor and cable away from heated surfaces. If either are damaged, they should be replaced with an approved Lumie product.

To power down or reset, remove the mains power adaptor from the wall socket. You should therefore ensure the socket is easily accessible.

Only to be used by adults and children over 7 years old. Children using the light, or playing near it, should be supervised.

Caution – do not attempt to service or modify any part of the unit. 

If the unit has been stored outside the recommended operating temperature range (5 °C to 35 °C / 41 °F to 95 °F), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.

The light source contained in this luminaire shall only be replaced by Lumie or an authorized Lumie distributor.

Not to be used by or on individuals who are unconscious, anaesthetized or incapacitated e.g. under the influence of alcohol.

CAUTION

Do not drop the unit on the floor, hit it hard or expose it to other shocks.

Place the unit on a stable, level and non-slippery surface.

Always unplug the unit before cleaning and allow it to cool down.

Not to be used with light dimmers, timers, motion detectors, voltage transformers or extension cords.

Intended for household use, including similar use in hotels. Do not use attachments not recommended by Lumie.

Disclaimer

It is not possible for Lumie to warn the consumer about every possible danger related to installation and use of the product. The consumer must use his or her own common sense and good judgment when installing and using the product.

Electromagnetic compatibility

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This Class B digital apparatus complies with Canadian ICES-005.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website:
www.lumie.com/pages/e-m-c.

Responsible Party Address

Violette Engineering Corporation
313 Park Avenue Suite 300
Falls Church, VA 22046 USA
info@violettecorp.com

Disposal

Do not throw out the mains power adaptor or the unit with normal domestic waste. Take them to an official recycling point to be disposed of responsibly.

Precautions

Research shows that light therapy is very safe. However, please check with your doctor first if any of the following apply:

- You have a major eye problem or have had eye surgery;
- You have been advised to avoid bright light or are taking medication that may cause photosensitivity;
- You have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can affect your mood.

A few people experience side effects but these are generally mild (including headaches, eye strain, nausea and hyperactivity) and can usually be solved simply by switching off the light. Over the next few days, gradually build up the time you spend in front of it to find a treatment time that works best for you.

Some people complain of difficulty sleeping and often this is down to using the light too close to bedtime.

If you have any problems while using Halo, or if you've been using the light for a week but haven't noticed any benefit, please contact Lumie or speak to your doctor.


IMPORTANT: the guidelines in this leaflet are not medical directions for treating any condition.

Introduction

Halo is a multifunctional light therapy solution, designed to look good and optimize your wellbeing am-to-pm, all year round.

The **main light** has adjustable brightness and color temperature, from bright white to help with daytime focus through to warmer tones for evening relaxation.

The **backlight** replicates the natural glow of dusk and dawn and allows you to create the perfect ambient lighting.

Main light in Day Mode

Day Mode utilizes both the warm-white and cool-white LEDs, allowing you to vary the color temperature and brightness to suit your needs. At maximum brightness, Halo delivers 10 000 lux at 20 cm / 8 in, simulating the natural light of the spring sun. Choose Day Mode to feel more energized, boost your mood, performance, alertness and productivity or to help treat the symptoms of winter blues.

Adjustable brightness

Adjustable color temperature: approx. 3000 K to 5000 K

Main light in Evening Mode

Evening Mode uses only warm-white LEDs with a lower color temperature, creating a more natural home setting after sundown. Even at full brightness, the reduced blue-light content of Halo in Evening Mode makes it ideal for reading, socializing or simply relaxing before bedtime.

Adjustable brightness

Fixed color temperature: approx. 3000 K

Backlight

The backlight provides further options for ambient lighting and can be used on its own or in addition to the main light. Warmer LEDs shine onto the copper-

colored reflector to produce a soft glow around the main light.

Fixed brightness

Fixed color temperature: approx. 2700 K

Please read these instructions carefully to get the most out of Halo.

Out of the box


Halo is supplied with two copper-colored legs.

- Slide the legs onto the two posts at the back of the unit.

The legs attach magnetically and allow for some movement to naturally support the unit as you set it down.

Take care not to lose the legs when lifting the assembled unit.

To remove the legs, simply pull firmly.

- Connect the mains power adaptor to the socket at the back of the unit.
- Position Halo on a flat surface.
- Plug in at the mains power.

The Lumie touch interface briefly illuminates whilst the touch sensor is calibrating.

Halo is now ready to use.

Note: If you are using Halo and the touch interface doesn't seem to be working properly, repeat the last steps and make sure you don't press the touch interface whilst the touch sensor is calibrating.

Lighting controls

The main light Day and Evening Modes are controlled by the touch interface:

Using Day Mode

- Touch .

This turns on the main light in Day Mode and also the backlight.

Day Mode offers variable color temperature and brightness to suit your needs throughout the day. At maximum brightness, the cool-white light is ideal for bright light therapy.

- Touch  to turn off both lights.

See also *Using the backlight*.

Using Evening Mode

- Touch .

This turns on the main light in Evening Mode and also

the backlight.

The reduced blue-light content of Evening Mode is warmer and ideal for use later in the day or close to bedtime. The maximum brightness is also reduced.

- Touch  to turn off both lights.

See also *Using the backlight*.

Switching between Evening and Day Mode

- Touch .

 is illuminated when in Evening Mode and you'll see that the light is more yellow.

Adjusting the brightness

- Touch and slide across **lumie**.

The letters illuminate from left to right to indicate increasing brightness.

In Day Mode this also adjusts the color temperature.

Using the backlight

- Touch  and  at the same time to turn off the backlight only.

The backlight is a fixed brightness light for ambient illumination. You can have the backlight on with or without the main light.

- Touch  and  at the same time to turn on the backlight only.

Day Mode for bright light therapy

Sunlight is a natural mood-booster, influencing body chemistry (e.g. melatonin and serotonin levels), making you feel more awake and positive. If you're not getting enough exposure to daylight, Halo is the next best thing and is available all year round.

- Make sure the main light is in Day Mode and at maximum brightness.
- Position Halo at arm's length (about 50 cm / 20 in) with the light shining towards your face.

Although you don't need to stare at it, the light must reach your eyes to have any effect. This means you can use Halo whilst working on your laptop, eating lunch, watching TV, reading, etc.

How often?

For wellbeing support

When you're indoors, use Halo in Day Mode to simulate sunlight. A bright light session works quickly and naturally to improve mood, energy, focus and performance whilst daily use can encourage a more regular sleep-wake pattern.

For winter blues

Look out for early signs of winter blues in Sept/Oct and if you start feeling more tired, anxious or irritable, start using Halo in Day Mode. As spring comes around, you can cut back on your bright light sessions. If symptoms return, simply use Halo again for a few days.

When?

If you're using Halo in Day Mode to help with the winter blues, bright light in the morning is most effective so try that first. For some people, early evening light exposure is beneficial.

I find it difficult to get going in the morning

Switch on Halo in Day Mode as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon.

I feel sleepy or fall asleep too early

Use Halo in Day Mode in the late afternoon/early evening. It's best to avoid bright light too close to bedtime – within about three hours – as you may find it difficult to get to sleep.

How long?

In Day Mode, with the light at full brightness, you're getting the maximum benefit so position Halo at arm's length (approx. 50 cm / 20 in) and start with 30 minutes a day (3000 lux).

You may find this is enough but if you don't notice much of an effect, increase your exposure each day until you find something that works for you.

If you're more comfortable with Halo positioned further away, then less light will reach your eyes and you'll need to use it for longer.

Similarly, if you're using Day Mode at reduced brightness, then less light will reach your eyes and you'll need to use it for longer.

There's no harm in using the light for longer than suggested. You can get your bright light therapy all in one go or use Halo on-and-off throughout the day for a cumulative effect.

Avoid using Day Mode at full brightness close to bedtime. Reduce the brightness or, better still, switch to Evening Mode with its less stimulating light tones.

If you find your eyes ache slightly, try using it further away, reduce the brightness or switch off for a while.

Whether you're using Halo for winter blues, or for general wellbeing, once you've found a pattern that's right for you try to use it around the same time every day.

Maintenance

Clean only with a dry or lightly moistened cloth.
Always unplug the unit first and let it cool if necessary.

Replacement mains power adaptors are available from Lumie.

LED lifetime

The LEDs in this product will typically last for many years, even with several hours use a day.

Warranty

This unit is covered under warranty by Lumie LLC against any manufacturing defect for 3 years from the date of purchase. If you need to return the product for servicing, use the original packaging (including internal packing) if possible. Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please contact Lumie for further advice.

For all warranty details, please refer to our Limited Warranty which can be found at <https://www.lumie.com/en-us/warranty>.

SN Serial number (see back of unit)

Technical description

Subject to change without notice.

Lumie Halo

Detachable accessory: mains power adaptor (supplied) 

Input: 100–240 V AC, 50/60 Hz, 0.85 A max.;

Output: 19 V $\approx$ 1.6 A (30.4 W)

Operating conditions:

Temperature: 5 °C to 35 °C / 41 °F to 95 °F

Relative humidity: 15% to 93% (non-condensing)

Atmospheric pressure: 70 kPa to 106 kPa

Storage/transport conditions:

Temperature: -25 °C to +70 °C / -13 °F to +158 °F

Relative humidity: 15% to 93% (non-condensing)

Atmospheric pressure: 25 kPa to 106 kPa

The weights, dimensions, capacities, performance ratings, characteristics and other data displayed are mentioned only as general information. They shall not bind Lumie.

Contact us

If you have any questions or problems regarding this product, or have suggestions that you feel could help us improve our products or service please email the Lumie Careline at info@lumie.com.

Distributed in the US and Canada by:
Lumie LLC, 4600 140th Avenue N., Suite 180,
Clearwater, FL 33762.

Online guide

You can view or download an online copy of the latest version of this guide on the appropriate product page of the Lumie website at <https://www.lumie.com/en-us/>.

Lumie retains the right to modify the guide as it deems appropriate.

The consumer is charged with the knowledge of the modifications set forth in the guide.

 Outside In (Cambridge) Limited, 3 The Links,
Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:
Outside In (Cambridge) Limited, incorporated in
England and Wales.
Registered Number: 2647359

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© Lumie 2025
Outside In (Cambridge) Limited, constituida en
Outside In (Cambridge) Ltd