lumie[®]

bodyclock[®]

shine 300

Safety

If the unit is damaged, contact Lumie before use.

Keep away from water and damp.

For indoor use only.

The unit gets warm when the light is on, so make sure that it is not covered.

Use only the mains power adaptor supplied.

Keep the mains power adaptor and cable away from heated surfaces. If either are damaged, they should be replaced with an approved Lumie product.

This appliance does not have a on/off switch. To power down or reset, remove the mains power adaptor from the wall socket. You should therefore ensure the socket is easily accessible.

Only to be used by adults and by children when supervised.

Do not attempt to service or modify any part of the unit.

If the unit has been stored outside the recommended operating temperature range (5 °C to 35 °C), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.

Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.

Electromagnetic compatibility

This product conforms to the electromagnetic compatibility (EMC) standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference. In rare cases where product function is lost due to EMC effects, disconnect the product from the power source for 5 seconds and then restart. If the problem persists, please contact Lumie.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Disposal

Do not throw out the mains power adaptor or the unit with normal domestic waste. Take them to an official recycling point to be disposed of responsibly.

	Serial number		Read instructions for use		Double-insulated (Class II electrical)
	European Declaration of Conformity: Medical Device Class I		Keep dry		DC (direct current) power supply
	Medical Device		For indoor use only		Storage/transport relative humidity
	Unique Device Identifier		Caution		Storage/transport temperature
	Manufacturer		Do not place in household waste. Dispose of at a recycling centre.		UK Conformity Assessed
	Authorised representative in the European Community		Member of the 'Green Dot' recycling scheme		
	Importer		Recyclable		

Introduction

Lumie Bodyclock wakes you naturally with a sunrise. Waking with light has been shown to improve mood, energy, productivity and the quality of sleep and awakening.¹

Wake-up lights can also help sufferers of SAD (Seasonal Affective Disorder) and winter blues cope with dark winter mornings.²

At bedtime, a fading sunset helps you naturally unwind by promoting the production of the sleep hormone, melatonin.³

Lumie invented the very first wake-up light in 1993. Use all year round to encourage a regular sleep-wake pattern.

Research shows that dawn simulation is very safe and without serious side-effects. If you have any concerns, please contact Lumie or speak to your doctor.

1 Acta Psychiatr Scand. 1993 Jul;88(1):67-71.

2 J Affect Disord. 2015 Jul 15;180:87-9.

3 Sleep Biol Rhythms. 2016 Jan;14(1):47-53.

Getting started

The ideal place for Lumie Bodyclock is on a bedside table near a mains socket.

- Connect the mains power adaptor to the socket in the base of the unit.
- Plug in the unit and switch on at the mains.

A  cursor will prompt you to enter the essential settings:

LANGUAGE

GENERAL		
LANGUAGE	ENG	
DAY	MON	
TIME	13:00	

- Use + or – to choose a language.
- Press **select** to confirm and continue.

DAY

GENERAL		
LANGUAGE	ENG	
DAY	MON	
TIME	13:00	

- Use + or – to set the day.
- Press **select** to confirm and continue.


TIME

GENERAL		
LANGUAGE	ENG	
DAY	MON	
TIME	10:56	

- Use + or – to set the hour.
- Press **select** to confirm.
- Do the same for the minutes.
- Press **select** to confirm and exit to the current time.

Bodyclock Shine 300 is now ready to use as a wake-up light, bedside light and radio.

Front panel

The four black buttons perform specific functions indicated by the icon next to them on the display. The bottom right icon indicates which audio option is in use.

To adjust the bedside light

Lumie Bodyclock can be used as a fully dimmable bedside light.

- Use $\wedge$ or $\vee$ to set the light level.
- Press and hold $\vee$ to turn off the light.

Timeout: If the bedside light is not turned off manually, the light will turn off automatically after 12 hours.

Setting the alarm time

Bodyclock Shine 300 has two types of alarm:

With DAILY you set one alarm time and turn it on manually when you go to bed.

With WEEKLY you can set the alarm time, or choose not to have an alarm, for each day individually. For example, you may want to wake up earlier on Monday and have no alarm at the weekend.

To choose DAILY or WEEKLY alarm type

- Press the settings button .

The display shows the first settings menu, ALARM:

ALARM TYPE

- Press **select** to edit.

The  cursor changes to  to show that you are editing the item.

- Use + or – to choose DAILY or WEEKLY.
- Press **select** to confirm and continue.

To set the DAILY alarm

If you selected DAILY use SET TIME to set an alarm time:

SET TIME

- Press **select** to edit.
- Use + or – to set the hour.
- Press **select** to confirm.
- Do the same for the minutes.
- Press the button next to  to exit from settings.

To set the WEEKLY alarm

If you selected WEEKLY you can choose for each day whether the alarm is OFF or ON, and if ON, the alarm time:

MON

- Press **select** to edit the day; for example, MON.
- Use + or - to choose OFF or ON and press **select** to confirm.

SET TIME

If you selected ON, an additional SET TIME item appears below the day.

- Use + or - to set the hour and press **select** to confirm.
- Do the same for the minutes.
- Repeat for each day.
- Press the button next to  to exit from settings.

After setting the alarm time(s), the next alarm is turned on and shown on the top line of the display:

Bodyclock Shine 300 is preset with a 30-minute sunrise and 30-minute sunset. To use it with these settings read the following sections.

To add an audio alarm or change any other settings see *Settings menus*, page 13, *Sunrise options*, page 15, and *Sunset options*, page 17.

To cancel an alarm

- Press 

OFF appears at the top of the display.

If you selected WEEKLY, when the cancelled alarm time has passed the next alarm time is shown.

Using the sunrise alarm

To wake with a sunrise alarm

If you set the alarm type to DAILY, you must activate the alarm before you go to sleep each time:

- Press the alarm button 

If you set the alarm type to WEEKLY, your alarm(s) will be turned on automatically.

The top line of the display shows the alarm time; for example:

In the morning

You'll wake with a gradual sunrise and any other sunrise options you have set. You can choose to have an internal sound or play the radio at the alarm time; see *Sunrise options*, page 15.

To snooze

- Tap the top of the lamp cover.

The sunrise, and audio alarm if set, will repeat after 9 minutes.

To deactivate the sunrise alarm

- Press any button.

This also turns off the audio alarm if you've selected one.

- Hold down  to turn off the light.

Timeout: If the light and audio alarm are not turned off manually, they will turn off automatically after 90 minutes.

Using the fading sunset

When you're ready to sleep, a fading sunset can help you to relax and unwind. You can choose to play an internal sound or the radio during the sunset; see *Sunset options*, page 17.

To start the sunset

- Use $\wedge$ or $\vee$ to set the light level.
- Press the sunset button .

The display shows  to confirm the sunset is activated.

The light will start to fade and any other sunset options you have set will be activated.

To stop the sunset

- Press the sunset button .


The display shows .

Display brightness

The display automatically adjusts to match the light level in your room when the room is dark and drops to low brightness when you're ready to sleep.

You can change this setting if you prefer the display to turn off completely when the room is dark.

- Use the switch in the base of the unit to set OFF or LOW:

To check the time when the display is off


- Tap the top of the lamp cover.

The display will change to low brightness for a few seconds.

Using the radio and internal sounds

Use + or – to set the volume, and < or > to tune the radio or select one of the internal sounds.

The bottom right icon shows which audio option is in use:

To select your audio option

- Press the audio button to choose FM radio **FM**, internal sounds  or no audio .
- Use + or – to set the volume.

The following sections give more information about the audio options.

To play internal sounds

Bodyclock Shine 300 has 15 sleep and wake sounds:

Beep, Tropical birds, Blackbird, Osprey, Goats, Guzheng, Cafe, Steam train, Thunderstorm, Waves, Crickets at night, Kitten purring, Tree frogs, Ping pong, White noise.

- Press the audio button until the icon next to it shows .
- Use < or > to choose a sleep/wake sound.

To use the FM radio

- Press the audio button until the icon next to it shows **FM**.
- Use < or > to tune manually.
- Press and hold < or > to search for the next station.

The station is shown on the bottom line of the display:

The signal strength is shown by the  icon. If necessary position the aerial for best reception.

To save a station as a preset

You can save up to five FM stations for instant listening:

- Use < or > to tune to the station you want to save.
- Press and hold **preset** until the display shows P1.
- Use < or > to choose the preset you want to use, then press **select**.

The display shows SAVED.

To listen to a saved station

- Press **preset** and use < or > to choose a saved station, then press **select**.

To turn off the audio

- Press the audio button so the icon next to it shows .

Settings menus

This section summarises the alarm, sunrise, sunset and other options available with Bodyclock Shine 300.

To access the settings menus

- Press the settings button .

The display shows the first settings menu, ALARM:

- Press < or > to step through the different menus.
- Press ^ or v to browse the current menu.

The currently highlighted item is shown by a  cursor to its right.

Longer menus will scroll to show you all the items.

To edit a setting

- Press **select**.

The  cursor changes to  to show that you are editing the item.

- Use + or – to edit the item.
- Press **select** to confirm the edit.

The  cursor will move to the next item on the menu.

To exit from the settings menus

- Press the button next to .

The following table summarises each menu and menu item:

Menu	Menu items	Options	Default
ALARM	ALARM TYPE	DAILY, WEEKLY	DAILY
	MON to SUN *	OFF, ON	ON
	SET TIME	HH:MM	07:00
SUNRISE	DURATION	OFF or 15 to 90 MIN	30 MIN
	LIGHT LEVEL	1 to 20	16
	ALARM SRC	Alarm source: OFF, RADIO or SOUND	OFF
	STATION *	P1 to P5	
	TRACK *	Filename	
	VOLUME *	1 to 20	10
SUNSET	DURATION	15 to 90 MIN	30 MIN
	NIGHTLIGHT	OFF or 1 to 10	OFF
	AUDIO SRC	OFF, MANUAL or SOUND	MANUAL
	TRACK *	Filename	
	VOLUME *	OFF, ON or FADE	OFF
GENERAL	LANGUAGE	ENG FRA ESP DEU NED SVE FIN NOR DAN ITA or POL	ENG
	DAY	MON to SUN	MON
	TIME		13:00

* Depending on previous settings.

To change the current time


Bodyclock Shine 300 will keep time for 30 minutes without power. All other settings are permanently saved.

If you need to set the time use the TIME setting in the GENERAL menu.

Sunrise options


- Press the settings button .
- Use < or > to choose the SUNRISE menu.
- Continue through the following sequence to set the sunrise options:

DURATION

- Press **select** to edit.
- Use + or - to set the duration of the sunrise: OFF or between 15 and 90 minutes.
- Press **select** to confirm and continue.


LIGHT LEVEL

Only available if you have set a DURATION.


- Press **select** to edit.
- Use + or - to set the final brightness of the sunrise from 1 to 20.
- Press **select** to confirm and continue.

ALARM SRC

- Press **select** to edit.
- Use + or - to set the audio alarm option for the end of the sunrise: OFF for no audio alarm, or RADIO, or SOUND to choose one of the wake sounds.
- Press **select** to confirm and continue.


STATION

Only available if you have set RADIO.

- Press **select** to edit.
 - Use + or - to set the radio preset for the end of the sunrise.
 - Press **select** to confirm and continue.
-

TRACK

Only available if you have set SOUND.

- Press **select** to edit.
- Use + or – to set the wake sound for the end of the sunrise.
- Press **select** to confirm and continue.

VOLUME

Only available if you have set RADIO or SOUND.

- Press **select** to edit.
- Use + or – to set the volume of the audio alarm from 1 to 20.
- Press **select** to confirm.
- Press the button next to  to exit from settings.

Sunset options

- Press the settings button .

- Use < or > to choose the SUNSET menu:

- Continue through the following sequence to set the sunset options:

DURATION

- Press **select** to edit.
- Use + or - to set the duration of the sunset: between 15 and 90 minutes.
- Press **select** to confirm and continue.

Note: The actual duration may be less than this if you start the sunset with the light already dimmed.

NIGHTLIGHT

- Press **select** to edit.
- Use + or - to set a nightlight: OFF for no nightlight (sunset fades to off) or a level from 1 to 10 (sunset fades to nightlight).
- Press **select** to confirm and continue.

AUDIO SRC

- Press **select** to edit.
- Use + or - to control any audio at the start of the sunset: OFF to cancel any audio, MANUAL to leave any audio playing or SOUND to choose one of the sleep sounds.
- Press **select** to confirm and continue.

TRACK

Only available if you have set SOUND.

- Press **select** to edit.
- Use + or - to set the sleep sound for the sunset.
- Press **select** to confirm and continue.

VOLUME

Only available if you have set MANUAL or SOUND.

- Press **select** to edit.
- Use + or – to control any audio at the end of the sunset: OFF to cancel any audio, ON to leave any audio playing or FADE to gradually decrease volume and turn off.
- Press **select** to confirm.
- Press the button next to  to exit from settings.

Maintenance

Clean only with a dry or lightly moistened cloth. Always unplug the unit first and let it cool if necessary.

Replacement mains power adaptors are available from Lumie.

To restore the default factory settings

- Turn off the power at the mains.
- Press and hold **select** and turn the power back on.
- Release the button when the software version number shows.

LED lifetime

The LEDs in this product will typically last for many years, even with several hours use each day. Although it's unlikely, it's possible for individual LEDs to fail. If this happens, contact Lumie; do not attempt to replace them.

Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you need to return the product for servicing, use the original packaging (including internal packing) if possible. Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please contact Lumie for further advice.

Note: Any serious incident relating to the unit should be reported to Lumie and to the competent authority in the country where you live.

Technical description

Subject to change without notice.

Bodyclock Shine 300

Detachable accessory: Mains power adaptor (supplied) 

Input: 100–240 Vac, 50/60 Hz, 1 A max.; Output: 15 V $\approx$ 1 A (15 W)

Operating conditions:

Temperature: 5 °C to 35 °C

Relative humidity: 15% to 93% (non-condensing)

Atmospheric pressure: 70 kPa to 106 kPa

Storage/transport conditions:

Temperature: -25 °C to +70 °C

Relative humidity: 15% to 93% (non-condensing)

Atmospheric pressure: 25 kPa to 106 kPa

European Declaration of Conformity: Medical Device Class I

UK Conformity Assessed

Medical Device

Contact us

If you have any questions or problems regarding this product, or have suggestions that you feel could help us improve our products or service please email info@lumie.com.

 Outside In (Cambridge) Limited, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

Registered Number: 2647359.

VAT No:

GB 880 9837 71 (UK)

NL 826301381B01 (EU)

 Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta.

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© Lumie 2022
Outside In (Cambridge) Limited, rekisteröity Englannissa ja Walesissa.
Outside In (Cambridge) Ltd